

Ordering Information

Item Code	Pore Size(μm)	Diameter(mm)	Sterile	Pre - filter	Quantity / pk
SLGF0445N	0.45	4	No	No	100pcs
SLGF04100N	1.0	4	No	No	100pcs
SLGF04160N	1.6	4	No	No	100pcs
SLGF1345N	0.45	13	No	No	100pcs
SLGF13100N	1.0	13	No	No	100pcs
SLGF13160N	1.6	13	No	No	100pcs
SLGF2545N	0.45	25	No	No	100pcs
SLGF25100N	1.0	25	No	No	100pcs
SLGF25160N	1.6	25	No	No	100pcs
SLGF3345N	0.45	33	No	No	50pcs
SLGF33100N	1.0	33	No	No	50pcs
SLGF33160N	1.6	33	No	No	50pcs
SLGF5045N	0.45	50	No	No	25pcs
SLGF50100N	1.0	50	No	No	25pcs
SLGF50160N	1.6	50	No	No	25pcs